college park
TECHNOLOGY for the HUMAN RACE

technical instructions

التعليمات الفنية • Teknisk vejledning • Technische Anleitung
Τεχνικές οδηγίες • Instrucciones técnicas • Teknisket ohjeet
Instructions techniques • הוראות טכניות • Istruzioni tecniche
Technische instructies • Tekniskie instruksjoner • Instrukcje techniczne
Instruções técnicas • Instruções técnicas • Технические инструкции
Technické pokyny • Tekniska anvisningar • Teknik Talimatlar • 技术说明

FIGURE 1



FIGURE 2

LANGUAGES

EN-English	4
AR-Arabic	7
BG-Bulgarian	10
CS-Čeština	13
DA-Dansk	16
DE-Deutsch	19
ES-Español	22
FI-Suomen Kieli	25
FR-Français	28
HR – Hrvatski	31
HU-Magyar	34
IT-Italiano	37
JA-Nihongo	40
KO-Kugo	43
NL-Nederlands	46
NO-Norsk	49
PL-Polski	52
PT-Português	55
PT-Português Brasil	58
RO-Română	61
RU-Russkiy	64
SK-Slovenčina	67
SI-Slovensko	70
SV-Svenska	73
TR-Türkçe	76
UK-Ukrayins'ka	79
ZH-Zhōngwén	82

Package Contents

- (1) Soleus Foot (1) Foot Shell (1) CPI Sock
(1) 4-Hole Fastener Kit (1) Cosmetic Attachment Plate (CAP)

Tools Recommended

- (1) 4 mm Hex Key
(1) Foot Horn

This diagram is to help familiarize you with the unique parts of the Soleus. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components (Figure 1)

- | | | | |
|-------------------------|-----------------------------|----------------------------|------------------------|
| A. Housing | B. Heel Spring Plate | C. Toe Spring Plate | D. Heel Springs |
| E. Mid Spring | F. Toe Spring | G. Mid Bumper | H. Heel Bumper |
| I. Flexible Belt | • CPI Sock (not shown) | • Foot Shell (not shown) | |

PRODUCT DESCRIPTION

This prosthetic foot device is constructed with a housing, four composite springs, two bumpers, and a flexible belt. The composite springs are secured to the housing and flexible belt with fasteners.

INTENDED USE

The Soleus is a prosthetic foot designed to replace one or more functions of the biologic human foot.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- None known

TECHNICAL SPECIFICATIONS

FOOT SIZE	WEIGHT LIMIT	BUILD HEIGHT	FOOT WEIGHT*
21-24 cm	220 lbs / 100 kg	6.4 in / 16.3 cm	667 g
25-26 cm	250 lbs / 113 kg	7.1 in / 18.0 cm	
27-30 cm	275 lbs / 125 kg	7.6 in / 19.3 cm	

*26cm foot w/shell

GAIT MATCHING® GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

WEIGHT LBS	0-100	101-120	121-140	141-160	161-180	181-200	201-220	221-250	251-275
WEIGHT KG	0-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125
SIZE CM	21-24								
Moderate Impact	11	12	33	34	54	55	55	N/A	N/A
High Impact	12	33	34	54	55	55	N/A	N/A	N/A
SIZE CM	25-26								
Moderate Impact	12	12	32	33	54	55	75	76	N/A
High Impact	12	32	33	54	55	75	76	76	N/A
SIZE CM	27-30								
Moderate Impact	11	11	12	32	33	34	54	55	55
High Impact	12	12	32	33	34	54	55	55	55

MOUNTING

Use only high quality proximal endoskeletal components. Use of mounting fasteners other than those provided will void the warranty.

ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the Foot Horn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

STATIC ALIGNMENT (Figure 2)

For optimal function, balance the patient's weight evenly between the heel and toe.

- The Soleus was designed with a 3/8" (10 mm) heel rise.
- The load line divides the foot at 1/3 heel lever and 2/3 toe lever.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE
Firmer Toe Response	Plantarflex the Soleus or move load line posterior
Softer Toe Response	Dorsiflex the Soleus or move load line anterior
Firmer Heel Response	Dorsiflex the Soleus or move load line anterior
Softer Heel Response	Plantarflex the Soleus or move load line posterior

MOUNTING FASTENER GUIDELINES

Fastener length will vary with adapter choice. M6 x 16mm fasteners are for use with rotatable adapters, while M6 x 14mm fasteners are for use with non-rotatable adapters.

FASTENER SELECTION	POTENTIAL RESULT
Thread length too short	Inadequate thread engagement
Thread length too long	Bottoming-out or binding

CAUTION

- Apply Loctite® 242® on all 4-hole fastener threads.
- Torque evenly to 12 N·m (8.8 ft-lbs).

Failure to follow these guidelines may result in one or all of the following: noise, movement, damage, and catastrophic failure.

EXPOSURE TO MOISTURE

To maintain the lifespan of the Enviroshell® and CPI sock, always dry foot thoroughly after exposure to moisture. Rinse and dry if exposed to salt water or chlorinated water.

WARNING

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

WARRANTY INSPECTION/ MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composites and Adapter
- CPI Sock
- Foot Shell

Warranty Inspection Schedule for Soleus: Annually.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc., and the competent authority of the Member State in which the prosthetist and/or patient is established.

*Serious incident' is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following; (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - LABEL

FOOT SIZE	WEIGHT LIMIT (KG)	LABEL TEXT
21-24 cm	100	ISO 10328-P5-100 kg
25-26 cm	113	ISO 10328-P6-113 kg
27-30 cm	125	ISO 10328-P6-125 kg

ISO 10328 - "P" - "m"kg*)

***) Body mass limit not to be exceeded!**

For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.

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©2026College Park Industries, Inc. All rights reserved.